Shaping New Perspectives on Pulmonary Function



Airwave Oscillometry Fast, Easy, Portable

THORASYS

tremoflo Benefits

Tidal Breathing No forced expiration.

K

Respiratory Mechanics Info related to small airways. ³



Fast & Easy Test in a few min.



Adult & Pediatric Normative values provided.





Easy Daily Calibration Less than 30sec.



Compact & Portable Weighs only 700 g.



Pre/Post testing Automated Workflows.



The Importance of Small Airways

The tremoflo C-100 provides information related to small airway function.

"The small airways are frequently involved early in the course of [asthma and COPD] diseases, with significant pathology demonstrable often before the onset of symptoms or changes in spirometry and imaging." McNulty and Usmani, ECRJ 2014¹

"Peripheral airway impairment may be clinically relevant at all levels of asthma severity and control." Galant et al., AAAI 2017²



tremoflo[°] Airwave Oscillometry Fundamentals-

tremoflo measurements are fast and easy. The patient just breathes quietly.

To assess respiratory function, the tremoflo adds a gentle oscillatory wave to the patient's regular, quiet breathing. A short measurement duration of only 20 seconds allows three repetitions within a couple of minutes, even in patients who have difficulty performing spirometry.

Oscillometry **Outcomes**

Summary of clinician interpretations per published studies

A pair of two curves calculated from the raw data reflects the mechanical properties of the respiratory system in characteristic patterns. Several key outcome parameters are then derived from these curves.

Characteristic Patterns

- The **Resistance** curve shown to reflect R central and peripheral airways.³
- The **Reactance** curve shown to reflect X elasticity and peripheral airways. 4

Key Outcomes

Resistance at 5 Hz ¹	R ₅	Normal	†
Resistance change, 5 to 20 Hz ¹	R ₅₋₂₀	Close to zero	Close to zero
Reactance area ¹	AX	Normal	Normal
Reactance difference insp-exp. ²	X _{5.in-ex}	Close to zero	Close to zero





Small





Actual tremoflo Data

A simple green-to-red gauge scale clearly shows whether patients fall within or outside normative values that are available for pediatrics and adults.



Expert Opinions

Oscillometry offers "an in-office tool that can assess peripheral airway impairment even in young children and may detect airway obstruction earlier than spirometry." *Galant et al., AAAI 2017*

"Monitoring small airway function by [Oscillometry] can be useful in identifying patients who are at risk for losing asthma control, and in assisting with clinical decisions and treatment plans." *Bickel et al., Chest 2014*

Oscillometry "allows for the discrimination of inspiratory and expiratory resistance and reactance ... to help discriminate between asthma and COPD." *McNulty & Usmani, ECRJ 2014*

Oscillometry "has been found to be useful in measuring response to bronchodilators, such as salbutamol and ipratropium, in patients with asthma and COPD." *Bickel et al., Chest 2014*

tremoflo C-100

Airwave Oscillometry System (AOS)
Measurement Principle Oscillometry (Forced Oscillation Technique FOT)

Medsarement i melple	Osemethetry (Foreca Osemation Teerinique, FO		
Oscillator Technology	Breathe-through Vibrating Mesh (Patented)		
Measurement Modes	AOS: Pseudo-random noise		
Measurement Duration	20 sec (typ.), user adjustable 3 repetitions (as per guideline)		
Patient Interface	Bacterial/viral filter with integrated mouthpiece		
Dimensions & Weight	19 x 13 x 14 cm, 0.7 kg (handheld only) 21 x 14 x 24 cm, 1.7 kg (handheld & cradle)		
Performance	Meets and exceeds ERJ 2003, 22: 1026-1041		
Covered Markets	Canada, United States, Israel		

More information at www.tremoflo.com

- ¹ Hirsh et al., AAAI 2011
- ² Vassiliou et al., EUR RESPIR J. 1996
- ³ Usmani et al., RESPIRATORY MEDICINE 2016
- ⁴ Al-Alwan et al., AJRCCM 2014
- ⁵ Galant et al., AAAI 2017
- ⁶ Johnson et al., THORAX 2007

The content of this document is believed to be accurate at the time of release. However, THORASYS and its affiliates offer no guarantees, expressed or implied, in case of typographic or other errors. All specifications are subject to change without notice.

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed practitioner.

THORASYS

Thoracic Medical Systems Inc.

Head Office

Montreal, Quebec, Canada +1 800 855 THORASYS European Office Teningen (Freiburg), Germany +49 7641 967 93 53

thorasys.com